



LUSTER

Dentovations introduces LUSTER, an at-home tooth whitening system that whitens teeth up to 6 shades in less than 1 hour with no tooth sensitivity. LUSTER's patented technology offers similar results to in-office treatments, but with less concentration of the oxidizing agent (peroxide) and less intensity of the light source. The patented resin formula in the Super Whitening Gel and the formula of the Accelerator Mouth Rinse, significantly enhances the efficacy of the whitening process.

DENTOVATIONS
800-420-8990 ext. 18095
<http://link.argifocus.com/OZH-95>
Reader Service 095



NUPRO WHITE GOLD

NUPRO White Gold in-office and take-home whiteners from DENTSPLY Professional were developed with clinician input, and they are available in both hydrogen peroxide and carbamide peroxide formulations. In-office results of the 36% hydrogen peroxide averaged 6 shades whiter in 45 minutes. The 15% hydrogen peroxide take-home whiteners require 15 minutes wear time each day for 7 days and, on average, whiten teeth 5 shades.

DENTSPLY PROFESSIONAL
800-420-8990 ext. 18096
<http://link.argifocus.com/OZH-96>
Reader Service 096



BRITESMILE

BriteSmile from Discus Dental offers a safe, effective way to whiten patient's teeth in just one office visit. The proprietary hydrogen peroxide gel is specially formulated to respond to the gentle blue wavelengths of the BriteSmile light system. To help maintain patient's brighter smile, BriteSmile also offers a range of post-care teeth whitening products.


DISCUS DENTAL
800-420-8990 ext. 18097
<http://link.argifocus.com/OZH-97>
Reader Service 097



NITEWHITE ACP AND DAYWHITE ACP

Only Discus Dental's NiteWhite ACP and DayWhite ACP take-home whitening systems contain amorphous calcium phosphate (ACP), potassium nitrate, and fluoride to provide excellent whitening while also rebuilding enamel, reducing fadeback, improving luster, and reducing sensitivity. The patented dual-barrel syringes provide the necessary stability to maintain shelf life without refrigeration. Kits can be personalized at no additional cost.


DISCUS DENTAL
800-420-8990 ext. 18098
<http://link.argifocus.com/OZH-98>
Reader Service 098




FENDERMATE®

The World's Fastest Matrix!


- Takes less than 5 seconds to apply...




Inserts as a wedge



Convex Contact Point Maker




Close Cervical Marginal Adaption avoids overhangs



Flexible wing separates the teeth and firmly seats the cervical margin

Brought to you by



www.jsdental.com

Manufactured by DIRECTA AB - Uplands Väsby, Sweden - www.directadental.com
 FenderMate® is a registrate trademark owned by Directa AB. Design and patent pending

JS Dental Mfg., Inc. - PO Box 904 - Ridgefield, CT 06877 - 203-438-8832 - 1-800-284-3368 - FAX 203-431-8485
 Reader Service 032